

SUCTION TECHNOLOGY

Ardo Suction Pumps Master / Senator

Strengths

- Easy to use, just a main switch and a vacuum regulator.
- High reliability thanks to the piston-cylinder-technology.
- Low noise level due to the Heidrive Motor.
Master: 43 dBA on max. vacuum
Senator: 42 dBA on max. vacuum
- "3guard" safety system with hydrophobic bacterial filter and overflow protection in the jar.
- Fast vacuum build up due to powerful engine construction.
- Service requested once a year, beside that maintenace free, replace the filter regularly.
- 5 years warranty

SWISS MADE 



Ardo Master Senator – Strengths / 12.2017 / V01